

Highlights of ELISA Kits from Biomatik

Last Update: 2024-10-03

Biomatik offers an extensive list of over 29,000 high-quality ELISA Kits. Most ELISA kits can be shipped within 3-11 business days for door-to-door delivery worldwide. You may [click here](#) for frequently asked questions about using our ELISA kits.

Highlights of ELISA Kits from Biomatik

- **High Sensitivity and Low Detection Range:** Biomatik's ELISA kits offer exceptional sensitivity with their lowest detection ranges in the picogram range (1-100 pg/mL, 3-200 pg/mL), making them ideal for detecting low-abundance targets, such as trace proteins and antibodies, with precision.
- **Versatile Target Detection:** The featured ELISA kits can accurately detect various targets, including antibodies, low-concentration proteins, and small molecules, even in micro-volume samples (15-25 μ L). With high sensitivity, the kits allow detection within a wide dynamic range, often completing assays in under an hour, ensuring quick and reliable results for diverse research needs.
- **Comprehensive Sample Compatibility:** Biomatik's ELISA kits are compatible with a broad range of sample types, including serum, plasma, tissue homogenates, cell lysates, cell culture supernatants, and other biological fluids, offering flexibility for research in multiple biological contexts.
- **Ease of Use and Simplified Protocols:** Many of the kits are designed with user-friendly protocols, featuring pre-diluted detection reagents and minimal handling steps. This makes them easy to use, even for researchers with limited ELISA experience. This allows for fast setup and increased productivity in the lab.
- **Enhanced Specifications for Efficiency:** Some of our kits feature optimized protocols with shortened experimental times and enhanced component stability, ensuring reproducibility and reducing handling time without compromising the quality of results.
- **Cited in Peer-Reviewed Research:** Biomatik's ELISA kits have been widely cited in scientific publications, demonstrating their reliability and trusted performance in the academic and scientific community. Researchers can confidently use these kits, knowing they have been validated in various published studies.



The table below highlights several popular targets currently on the market. We encourage you to explore additional targets of interest on our website.

No.	Cat#	Product Name
1	EKU11047	Lipopolysaccharide (LPS) ELISA Kit
2	EKU03971	Estradiol (E2) ELISA Kit
3	EKU06215	Human Neurofilament, Light Polypeptide (NEFL) ELISA Kit
4	EKU06773	Human Procollagen I N-Terminal Propeptide (PINP) ELISA Kit
5	EKU07605	Testosterone (Testo) ELISA Kit
6	EKU07900	Human Cardiac Troponin I (cTnI) ELISA Kit
7	EKU05881	Human Matrix Metalloproteinase 9 (MMP9) ELISA Kit
8	EKU03476	Cortisol (Cor) ELISA Kit
9	EKU03779	Dopamine (DA) ELISA Kit
10	EKU05783	Malondialdehyde (MDA) ELISA Kit

Restrictions: All Products are for Research Use Only