

anti- RSK2 antibody

Product Information

Catalog No.:	CAF50308
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	ribosomal protein S6 kinase, 90kDa, polypeptide 3
Synonyms:	None
Observed MW:	84 kDa
Uniprot ID :	P51812

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human kidney tissue slide using CAF50308(RPS6KA3 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with CAF50308(RPS6KA3 antibody) at dilution of 1:1000