

anti- 5-HT antibody

Product Information

Catalog No.:	CAF50403
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen: 5-hydroxytryptamine

Application

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:50-1:500

Image:



Immunohistochemistry of paraffin-embedded rat brain tissue slide using CAF50403 (5-HT antibody) at dilution of 1:200

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.